



Analysis of TO-15/TO-17 air toxics in urban air using TD-GC/TOF MS and automated compound identification software


NEMC, Washington D.C. August 6th 2012



Nicola Watson
Environmental Specialist
Markes International
nwatson@markes.com



Thermal desorption – One versatile technique for all vapour-phase air monitoring applications



Application: 'Air toxics' in canisters – US EPA Method TO-15


- 1 Propylene
- 2 Dichlorodifluoromethane
- 3 1,2-Dichlorotetrafluoroethane
- 4 Methyl chloride
- 5 1,2-Dichloroethane
- 6 1,3-Butadiene
- 7 Vinyl chloride
- 8 Methyl bromide (bromomethane)
- 9 Chloroethane
- 10 Trichlorotrifluoroethane (Freon® 113)
- 11 Ethanol
- 12 1,2-Dichloroethylene
- 13 1,1,2-Trichlorotrifluoroethane
- 14 Acetone
- 15 Carbon disulfide
- 16 Isopropyl alcohol
- 17 Methylene chloride
- 18 Tert-butyl methyl ether
- 19 *n*-Hexane
- 20 1,1-Dichloroethane
- 21 Vinyl acetate

- 22 Cis-1,2-Dichloroethylene
- 23 Methyl ethyl ketone
- 24 Ethyl acetate
- 25 Tetrahydrofuran
- 26 Chloroform
- 27 1,1,1-Trichloroethane
- 28 Cyclohexane
- 29 Carbon tetrachloride
- 30 Benzene
- 31 *n*-Heptane
- 32 Trichloroethylene
- 33 1,2-Dichloropropane
- 34 1,4-Dioxane
- 35 Bromodichloromethane
- 36 Trans-1,3-dichloropropene
- 37 Methyl isobutyl ketone
- 38 Toluene
- 39 Cis-1,3-Dichloropropene
- 40 Trans-1,2-Dichloroethylene
- 41 1,1,2-Trichloroethane
- 42 Tetrachloroethylene

- 43 Methyl *n*-butyl ketone
- 44 Dibromochloromethane
- 45 1,2-Dibromoethane
- 46 Chlorobenzene
- 47 Xylene
- 48 Xylene
- 49 Xylene
- 50 Styrene
- 51 Tribromomethane
- 52 1,1,2,2-Tetrachloroethane
- 53 1,2,4-Trimethylbenzene
- 54 1,3,5-Trimethylbenzene
- 55 1-Ethyl-4-methyl benzene
- 56 Ethylbenzene
- 57 1,2-Dichlorobenzene
- 58 1,3-Dichlorobenzene
- 59 Chloromethylbenzene (alpha)
- 60 1,4-Dichlorobenzene
- 61 1,2,4-Trichlorobenzene
- 62 Hexachloro-1,3-butadiene




1 L of a 1 ppb air
toxics mix
analysed splitless
and cryogen-free
using TD-GC/MS
scan



Application: 'Air toxics' using sorbent tubes – US EPA Method TO-17

1	Propylene	22	Cis-1,2-Dichloroethylene	43	Methyl <i>n</i> -butyl ketone
2	Dichlorodifluoromethane	23	Methyl ethyl ketone	44	Dibromochloromethane
3	1,2-Dichlorotetrafluoroethane	24	Ethyl acetate	45	1,2-Dibromoethane
4	Methyl chloride	25	Tetrahydrofuran	46	Chlorobenzene
5	1,2-Dichloroethane	26	Chloroform	47	Xylene
6	1,3-Butadiene	27	1,1,1-Trichloroethane	48	Xylene
7	Vinyl chloride	28	Cyclohexane	49	Xylene
8	Methyl bromide (bromomethane)	29	Carbon tetrachloride	50	Styrene
9	Chloroethane	30	Benzene	51	Tribromomethane
10	Trichlorotrifluoroethane (Freon® 113)	31	<i>n</i> -Heptane	52	1,1,2,2-Tetrachloroethane
11	Ethanol	32	Trichloroethylene	53	1,2,4-Trimethylbenzene
12	1,2-Dichloroethylene	33	1,2-Dichloropropane	54	1,3,5-Trimethylbenzene
13	1,1,2-Trichlorotrifluoroethane	34	1,4-Dioxane	55	1-Ethyl-4-methyl benzene
14	Acetone	35	Bromodichloromethane	56	Ethylbenzene
15	Carbon disulfide	36	Trans-1,3-dichloropropene	57	1,2-Dichlorobenzene
16	Isopropyl alcohol	37	Methyl isobutyl ketone	58	1,3-Dichlorobenzene
17	Methylene chloride	38	Toluene	59	Chloromethylbenzene (alpha)
18	Tert-butyl methyl ether	39	Cis-1,3-Dichloropropene	60	1,4-Dichlorobenzene
19	<i>n</i> -Hexane	40	Trans-1,2-Dichloroethylene	61	1,2,4-Trichlorobenzene
20	1,1-Dichloroethane	41	1,1,2-Trichloroethane	62	Hexachloro-1,3-butadiene
21	Vinyl acetate	42	Tetrachloroethylene		

Splitless desorption of
'Air toxics' tube loaded
with 1 L of 1 ppb std
GC/MS




Can canisters do everything?



- ✓ Great for C₂ to C₁₂ compounds
- ✓ Suitable for rapid transfer (not storage) of ultra-volatile reactive compounds such as H₂S
- ✓ Ideal for simple grab-sampling

- ✗ NOT suitable for compounds with volatility less than C_{10/12}
- ✗ NOT suitable for high-concentration samples
- ✗ Time-weighted average sampling is NOT easy with a canister


Application: Soil gas




- Profiles of soil gas contaminated with kerosene obtained using:
- (a) Canister sampling and TO-15 analysis (blue)
 - (b) Sorbent tube sampling with TO-17 analysis (red)

Active (pumped) sampling

- TO17 specifies “The monitoring procedure involves pulling a volume of air through a sorbent packing to collect VOCs”
- Flow rate typically 20 – 100 ml/min
- Volume taken is typically in the range 5 ml – 100 L, depending on expected concentration (typically two samples 1 L and 4 L)
- Much faster technique compared to diffusive sampling
- Important do not exceed breakthrough volume for a compound on a given sorbent




Tube-based thermal desorption – An overview of the process




Tube-based thermal desorption

Sample passes onto the sorbent




Compounds of interest are adsorbed on the sorbent surface

Tube-based thermal desorption



Lighter gases such as nitrogen pass through the sorbent

Breakthrough



Affected by:


- | | |
|-----------------|--|
| Temperature | lower temp. = stronger interaction |
| Sample volume | lower volume = less risk of breakthrough |
| Mass or sorbent | more sorbent = more surface area |
| Type of sorbent | stronger sorbent = stronger interaction |

Air monitoring – Pumped

Sorbent selection for both tubes and focusing trap are very important


Semi-volatile compounds – Weak sorbent

Helps prevent retention of unwanted compounds



Very volatile compounds – Strong sorbent

Prevents **breakthrough** of light compounds




Common sorbents

Sorbent name	Volatility range
Quartz wool / silica beads	C ₃₀ – C ₄₀
Tenax TA	C ₇ – C ₃₀
Carbograph 2TD	C ₈ – C ₂₀
Carbograph 1TD	C _{5/6} – C ₁₄
Carbograph 5TD	C _{3/4} – C _{6/7}
SulfiCarb	C ₃ – C ₈
Carboxen 1003	C ₂ – C ₅
Carbosieve SIII	C ₂ – C ₅

Water retention




Application examples – High/Low concentration




50 mL sample of diesel exhaust

500 mL sample of 4 ppb ozone precursor standard



Sample security using sample re-collection


Stage 1: Primary (tube) desorption with optional (inlet) split



- Patented heated valve is inert and low volume: Allows quantitative recovery of high and low volatility and reactive compounds
- The heated valve isolates the TD system allowing method compliance: leak testing, backflush trap desorption, purge to vent, overlap mode, etc.

Sample security using sample re-collection


Stage 2: Secondary (trap) desorption with optional (outlet) split



- Repeat analysis of re-collected samples makes it easy to validate analyte recovery through the TD flow path
- A change to the overall VOC profile indicates any bias

Re-analysis of low-concentration sample

The 4 ppb standard was re-collected for re-analysis using SIM detection conditions



BenchTOF-dx: Detector enhancements for air monitoring





What does BenchTOF-dx offer?

- **Spectral accuracy** cannot be compromised
- **Sensitivity** is **KING**
- **Speed** can be leveraged for deconvolution
- **Selectivity** – enhanced mass resolution should mainly be used to limit the matrix in VOC work (high res has limited advantages)
- **Stability** is key to productivity

Quadrupole comparison


500 mL sample of 4 ppb ozone precursor standard



Detection method	S/N
Full scan (Quad)	15:1
SIM (Quad)	200:1
BenchTOF-dx	1500:1


Quadrupole comparison

200 mL sample of ambient rural air



Quadrupole comparison

200 mL sample of ambient rural air




Carbon tetrachloride

Atmospheric concentration
~100 ppt (~ 85 pg on column)

Detection method	S/N
Full scan (Quad)	ND
SIM (Quad)	100: 1
BenchTOF-dx	700: 1

10 mL of ambient semi-rural air



Total ion chromatogram showing splitless analysis of only 10 mL of semi-rural air using TD-GC-TOF MS.

Inset: Extracted-ion chromatogram for a characteristic fragment ion of Freon® 113 (present in the atmosphere at ca. 80 ppt).

How can I use a large sensitivity boost in air monitoring applications?


- Trace-level work for unknowns and targets combined at lower MDLs
- Smaller sample sizes but same MDLs
- Higher splits, cleaner system but same MDLs
- However you want to!

An investment in BenchTOF-dx provides a sensitivity boost!

Provides productivity too!


Quad MS – 1 mL injection of 1 ppm standard
(62 component) with no split

Equivalent to ~1000 ppt (1 ppb) on column



BenchTOF-dx – 1 mL injection of 1 ppm standard
(65 component) with 100:1 split


Equivalent to ~10 ppt (0.01 ppb) on column



...without compromising sensitivity

BenchTOF-dx – 1 mL injection of 1 ppm standard (65 component) with 292:1 split

Equivalent to ~3 ppt (0.003 ppb) on-column




RMS signal-to-noise
ranges from
15:1 to 1350:1
(non DBC)

LODs (assuming 3:1
minimum S/N)
0.01 ppt (10 ppq) to
0.6 ppt (600 ppq)

BenchTOF-dx – at least 100 times better S/N than a quadrupole in
full scan mode

...whilst maintaining linearity



Challenges associated with identifying target compounds in complex GCMS profiles

The problem...

- Identifying the presence of known toxic chemicals in complex GCMS TIC profiles is very challenging

The way forward....


- Technical capacity to interpret spectra is declining or absent, majority of analysts rely on commercial libraries (NIST, Wiley etc)
- Development of data mining software which looks more intensively at the GCMS data
- Implementation of spectral deconvolution to deconvolve complex analyte/matrix spectra and matrix effects, producing mixed spectra
- The application of multivariate data analysis algorithms provides an orthogonal dimension of analysis. Conventional library search techniques can result in mis identification with low quality matches and limited confidence in the result



MARKES
international

Data Analysis using TargetView software

Target compound	CAS no.	Retention time (mins)	Expected retention time (mins)	delta RT (seconds)	Retention index library	Matching coefficient	Peak sum (TIC)	Peak sum (extr. ion)	Extracted ion
2-Butanone	78-93-3	10.504	-	-	-	0.824	737923	-	-
Acetone	67-64-1	7.084	-	-	-	0.916	5385525	-	-
Benzene	71-43-2	12.044	-	-	-	0.949	2806190	-	-
Benzene, 1,3-dichloro-	541-73-1	21.635	-	-	-	0.919	650369	-	-
Benzene, 1,3-dimethyl-	108-38-3	17.972	-	-	-	0.947	2839083	-	-
Benzene, 1-ethyl-4-methyl-	622-96-8	20.073	-	-	-	0.817	120012	-	-
Dichlorodifluoromethane	75-71-8	3.493	-	-	-	0.933	2271907	-	-
Ethane, 1,1,2-trichloro-1,2,2-trifluoro-	76-13-1	6.972	-	-	-	0.913	810621	-	-
Ethanol	64-17-5	6.164	-	-	-	0.892	1065482	-	-
Ethylbenzene	100-41-4	17.767	-	-	-	0.951	2826414	-	-
Heptane	142-82-5	12.442	-	-	-	0.940	281461	-	-
Hexane	110-54-3	9.114	-	-	-	0.949	808826	-	-
Isopropyl Alcohol	67-63-0	7.364	-	-	-	0.864	1471031	-	-
Methylene Chloride	75-09-2	8.082	-	-	-	0.865	378759	-	-
Styrene	100-42-5	18.682	-	-	-	0.943	9980339	-	-
Tetrachloroethylene	127-18-4	16.102	-	-	-	0.970	155316	-	-
Toluene	108-88-3	15.271	-	-	-	0.948	4746104	-	-
Trichloroethylene	79-01-6	13.129	-	-	-	0.918	191774	-	-
Trichloromethane	67-66-3	7.760	-	-	-	0.814	52008	-	-



Chromatogram of semi-rural air sample, before and after DBC processing. The inset shows TargetView report of those TO-15/17 compounds positively identified in the sample

Dynamic Background Compensation (DBC)


- **Dynamic Background Compensation (DBC) is a sophisticated software algorithm designed to distinguish between chromatographic peaks and slower GC-background / baseline anomalies.**

The main advantages of DBC include:

- **Improved spectral purity**
 - Enhanced identification of trace target analytes and unknowns.
- **Selective elimination of interfering ions resulting in**
 - Flat chromatographic baseline
 - Enhanced integration
- **Increased sensitivity**
 - Reduced noise enhances S/N


DBC in action –

Baseline noise suppression




Original baseline optimized, minimal background (TIC) is a

Deconvolution of peaks



Deconvolved spectra for both known and unknown compounds can be cross searched against NIST databases

The chemometric visualisation of data (PCA)



PCA identification of benzene

The close proximity of 1 spectrum indicates a single positive match

- Deconvolution identifies a number of target and unknown compounds in the sample
- The deconvolved spectrum for all compounds (yellow circles) is projected onto a 2 dimensional plane.
- The red dot equates to the target compound spectrum
- Proximity of a yellow circle to the target indicates a positive hit

Simplicity of use -

Generating a TargetView report is just a 2 step process

- 1) Select a data file
- 2) Import the data and print a report

Import Chemstation File

Chromatogram: Pesticides_Orangeoil_40_dbc.d

Select: Browse for Folders (*.d)

Peak view

Report File

Library: Pesticides.MSP

Chromatogram: Pesticides_Orangeoil_40_dbc

Sort by: Retention time

Show hits only

Target compound	CAS no.	Retention time (mins)	Expected retention time (mins)	delta RT (seconds)	Retention index library	Matching coefficient	Peak sum (TIC)	Peak sum (extr. ion)	Extracted ion
o-Hydroxybiphenyl	90-43-7	3.418	-	-	-	0.923	67943205	-	-
Diazinone	333-41-5	5.400	-	-	-	0.888	75340104	-	-
Chlorpyrifos - methyl	5598-13-0	6.168	-	-	-	0.834	73567104	-	-
Methyl parathion	298-00-0	6.284	-	-	-	0.904	52940385	-	-
Pirimiphos methyl	29232-93-7	6.753	-	-	-	0.857	64755360	-	-
Chlorpyrifos	2921-88-2	7.067	-	-	-	0.884	65482245	-	-
Methodathion	950-37-8	8.124	-	-	-	0.888	66934715	-	-
Ethion	563-12-2	9.180	-	-	-	0.891	70387200	-	-
Chlorobenzilate	510-15-6	9.681	-	-	-	0.895	4771600	-	-
Propargite	2312-35-8	9.766	-	-	-	0.786	61335930	-	-
Bromopropylate	18181-80-1	10.092	-	-	-	0.795	46128300	-	-
Benzene, 1,2,4-trichloro-5-[(4-chlorophenyl)sulfonyl]-	116-29-0	10.329	-	-	-	0.893	68813394	-	-

PRINT EXIT

Import options:

Minimum signal of most abundant ion [counts]: 1000

Deconvolution:

Automatically process target library:

Target library: Pesticides.MSP

TIC: Ion trace:

Select Cancel

Import Chromatogram Cancel

Summary

- A **combination** of canisters and sorbent tubes provides a comprehensive evaluation of an application, e.g. ambient air monitoring
- **BenchTOF-dx** provides a sensitivity and productivity boost that can be utilised several ways while providing method tunes and NIST-compliant spectra
- Challenges associated with identifying target compounds in complex GCMS profiles can now be overcome by employing simple data analysis packages such as **TargetView**.

Any Questions?

Nicola Watson
Environmental Specialist
Markes International
nwatson@markes.com

