



HUMAN HEALTH | ENVIRONMENTAL HEALTH




## Trace Metals Analysis of Flue Gas Desulfurization Wastewaters by ICP-MS

Stan Smith and Ewa Pruszkowski, Ph.D.,  
PerkinElmer, Inc.

NEMC 2012

- ▶ In 2009 US EPA completed study of wastewater discharges from the steam electric power generating industry
- ▶ Focus was on wastewater from flue gas desulfurization (FGD) systems as these comprised a significant portion of the plant pollutant discharges
- ▶ Rulemaking process was begun to address pollutants and waste streams not already covered in 40 CFR Part 423 (Steam Electric Power Guidelines)
- ▶ Proposed revisions to the effluent guidelines are expected November of 2012 and final rule expected April 2014

Typical FGD treatment system using wet scrubber to remove SO<sub>2</sub> emissions from the flue gas



- ▶ ICP-MS is currently the preferred technique for FGD due to its sensitivity and wide availability
- ▶ Current methods such as 200.8, 6020, and 1638 lack detailed treatment for interferences associated with FGD
- ▶ EPA has proposed a new draft test method: "Standard Operating Procedure: Inductively Coupled Plasma-Mass Spectrometry for Trace Element Analysis in Flue Gas Desulfurization Wastewaters"
- ▶ Just what does FGD wastewater look like?

- ▶ Constituents of FGD wastewater come from the limestone used to make the slurry, the coal burned, and the water supply
- ▶ FGD wastewaters vary widely depending on above parameters as well as the efficiency of the FGD treatment system

Parameter	Range (mg/L)
Total Dissolved Solids	4000 - 68,000
Chloride	1000 - 35,000
Sulfate	1500 - 8000
Calcium	750 - 4000
Magnesium	1100 - 4800
Sodium	670 - 4800
Boron	10 - 800
Total Organic Carbon	5 - 1100

- ▶ Evaluate Draft SOP
- ▶ Efficacy of Digestion Procedure (1638) and Possible Contamination from Open Vessel Digestion vs Closed Vessel Digestion
- ▶ Necessity of Dilution (1:10)
  - If on-line dilution, what conc. for Cal, QC, and Synth FGD?  
What about Ag?
- ▶ Apply both Collision Mode with KED as well as Reaction Mode using Perkin Elmer NexION model 300D ICP-MS with Universal Cell Technology
- ▶ Discreet sample intro from ESI; SC2 DX Autosampler and FAST valve

# Instrument Conditions

Component/Parameter	Type/Value/Mode
Nebulizer	Concentric PFA
Spray Chamber	Peltier-cooled baffled quartz cyclonic
Interface Cone Material	Nickel
Plasma Gas Flow	16.0 L/min
Auxiliary Gas Flow	1.2 L/min
Nebulizer Gas Flow	0.98 L/min
Sample Uptake Rate	270 $\mu$ L/min
RF Power	1600 W
Analytes	Al, As, Cd, Cr, Cu, Pb, Mn, Ni, Se, Ag, Tl, V and Zn
Internal Standards	Sc, Ge, In, (added on line)
Isotopes Monitored (details to follow)	22 (analytes and internal standards)
Modes of Operation (details to follow)	1. Standard 2. Collision/KED (He gas) 3. Reaction (NH <sub>3</sub> gas)
Replicates per Sample	3
Measurement time (3 reps)	1 min 30 sec
Analysis time (sample to sample)	2 min 30 sec

# Typical Analytical Sequence

Turn on instrument	Analyzed once per sequence
Warm-up instrument	
Tune instrument	
Perform mass calibration	
Perform resolution check	
Validate tuning criteria	
Calibration blank	
Calibration standards	
Initial Calibration Verification (ICV)	
Initial Calibration Blank (ICB)	
Reporting Limit (RL) verification standard	
Synthetic FGD wastewater matrix	
Spiked synthetic FGD matrix	
Continuing Calibration Blank (CCB) (carryover check)	
Continuing Calibration Verification (CCV)	
Continuing Calibration Blank (CCB)	Must be repeated
10 Samples, including laboratory control sample (LCS) plus one matrix spike (MS) and matrix spike duplicate (MSD)	
CCV	
CCB	



- ▶ Seemingly simple mixture
- ▶ Purchase custom solution from commercial supplier
- ▶ Prepare from neat reagents and stock elemental solutions
- ▶ Beware contamination in stock solutions

<b>Matrix</b>	<b>Concentration (mg/L)</b>
Chloride	5000
Calcium	2000
Magnesium	1000
Sulfate	2000
Sodium	1000
Butanol	2000

# Isotopes, Analysis Mode, IDL, and MDL

<b>Analyte/Mass (amu)</b>	<b>Mode</b>	<b>Internal Standard</b>	<b>* DL (µg/L)</b>	<b>MDL (µg/L)</b>
<b>Al - 27</b>	Standard	Sc	0.040	0.54
<b>V - 51</b>	Reaction	In	0.001	0.02
<b>Cr - 52</b>	Reaction	In	0.005	0.09
<b>Mn - 55</b>	Reaction	In	0.007	0.27
<b>Ni - 60</b>	Collision	Ge	0.005	0.28
<b>Cu - 63</b>	Collision	Ge	0.011	0.30
<b>Zn - 66</b>	Collision	Ge	0.065	1.20
<b>As - 75</b>	Collision	Ge	0.019	0.30
<b>Se - 78</b>	Collision	Ge	0.190	2.20
<b>Ag - 107</b>	Standard	In	0.001	0.03
<b>Cd - 111</b>	Collision	Ge	0.007	0.10
<b>Tl - 205</b>	Standard	In	0.001	0.01
<b>Pb - 208</b>	Standard	In	0.002	0.20

# Synthetic FGD Results

- ▶ 3 of several synthetic FGDs; 2 commercial and one in-house


<b>Analyte/Mass (amu)</b>	<b>Synthetic FGD #1 (µg/L)</b>	<b>Synthetic FGD #2 (µg/L)</b>	<b>Synthetic FGD #3 (µg/L)</b>
<b>Al - 27</b>	36.4	5.33	7.16
<b>V - 51</b>	5.11	4.31	0.19
<b>Cr - 52</b>	2.97	0.76	1.38
<b>Mn - 55</b>	37.3	25.4	7.02
<b>Ni - 60</b>	4.36	7.44	1.42
<b>Cu - 63</b>	2.35	1.78	3.03
<b>Zn - 66</b>	5.84	27.6	9.26
<b>As - 75</b>	1.11	2.88	1.09
<b>Se - 78</b>	3.10	4.21	3.04
<b>Ag - 107</b>	0.22	2.33	0.31
<b>Cd - 111</b>	3.25	3.92	0.11
<b>Tl - 205</b>	0.19	0.94	0.05
<b>Pb - 208</b>	1.36	9.23	5.37

# Spiked Synthetic FGD Results

<b>Analyte/Mass (amu)</b>	<b>Synthetic FGD #3 (µg/L)</b>	<b>Spike Value (µg/L)</b>	<b>Spike Recovery (%)</b>
<b>Al - 27</b>	7.16	40	111
<b>V - 51</b>	0.19	40	114
<b>Cr - 52</b>	1.38	40	107
<b>Mn - 55</b>	7.02	40	112
<b>Ni - 60</b>	1.42	40	100
<b>Cu - 63</b>	3.03	40	93.0
<b>Zn - 66</b>	9.26	40	94.0
<b>As - 75</b>	1.09	40	110
<b>Se - 78</b>	3.04	40	108
<b>Ag - 107</b>	0.31	40	90.0
<b>Cd - 111</b>	0.11	40	107
<b>Tl - 205</b>	0.05	40	103
<b>Pb - 208</b>	5.37	40	94.2

# Continuing Calibration Verification Stability


- ▶ Plot of CCV results for 110 sample sequence



# Internal Standard Performance

- ▶ Plot of internal standards for entire run

“Clean” matrix of CCV and CCB



Suppression from FGD matrix

# Real World FGD Matrix Spikes

<b>Analyte/ Mass (amu)</b>	<b>Sample #1 (µg/L)</b>	<b>Sample #1 + Spike (µg/L)</b>	<b>Sample #1 + Spike Dup</b>	<b>RPD</b>	<b>Spike Recovery (%)</b>
<b>Al - 27</b>	18.5	58.5	59.1	1.0	99.9
<b>V - 51</b>	1.46	47.1	48.6	3.1	114
<b>Cr - 52</b>	1.59	45.1	45.7	1.3	109
<b>Mn - 55</b>	10.1	53.7	54.6	1.6	109
<b>Ni - 60</b>	0.33	39.8	39.2	1.5	98.6
<b>Cu - 63</b>	0.11	38.4	38.2	0.5	95.6
<b>Zn - 66</b>	1.14	41.4	42.3	2.1	101
<b>As - 75</b>	0.50	40.9	40.9	0.1	101
<b>Se - 78</b>	47.8	100.8	99.4	1.4	132
<b>Ag - 107</b>	<0.003	34.9	33.7	3.3	87.1
<b>Cd - 111</b>	<0.01	40.6	40.8	0.5	101
<b>Tl - 205</b>	1.10	38.6	37.5	2.7	93.7
<b>Pb - 208</b>	0.007	34.8	33.7	3.1	86.9

# Real World FGD Matrix Spikes

<b>Analyte/ Mass (amu)</b>	<b>Sample #2 (µg/L)</b>	<b>Sample #2 + Spike (µg/L)</b>	<b>Sample #2 + Spike Dup</b>	<b>RPD</b>	<b>Spike Recovery (%)</b>
<b>Al - 27</b>	16.6	55.9	56.6	1.3	98.3
<b>V - 51</b>	0.19	46.0	44.8	2.8	115
<b>Cr - 52</b>	0.19	43.6	43.2	0.8	108
<b>Mn - 55</b>	1.81	44.7	45.0	0.7	107
<b>Ni - 60</b>	4.96	45.6	45.0	1.2	101
<b>Cu - 63</b>	2.24	39.5	38.4	2.7	93.1
<b>Zn - 66</b>	1.74	40.4	38.8	3.9	96.6
<b>As - 75</b>	0.40	44.8	44.6	0.5	111
<b>Se - 78</b>	36.5	84.2	81.1	3.7	119
<b>Ag - 107</b>	<0.003	35.6	37.2	4.5	89.0
<b>Cd - 111</b>	<0.01	42.0	40.4	3.9	105
<b>Tl - 205</b>	0.03	42.1	43.6	3.6	105
<b>Pb - 208</b>	0.09	39.3	41.2	4.6	98.1



# Real World FGD Matrix Spikes

<b>Analyte/ Mass (amu)</b>	<b>Sample #3 (µg/L)</b>	<b>Sample #3 + Spike (µg/L)</b>	<b>Sample #3 + Spike Dup</b>	<b>RPD</b>	<b>Spike Recovery (%)</b>
<b>Al - 27</b>	4.82	40.2	40.9	1.7	88.5
<b>V - 51</b>	0.04	44.3	44.1	0.4	111
<b>Cr - 52</b>	0.09	42.3	42.7	1.0	105
<b>Mn - 55</b>	2191	2194	2227	1.5	N/A
<b>Ni - 60</b>	0.74	40.3	39.3	2.3	98.8
<b>Cu - 63</b>	0.04	35.6	35.3	0.7	88.9
<b>Zn - 66</b>	0.42	36.3	37.1	2.1	89.7
<b>As - 75</b>	0.60	43.8	43.5	0.7	108
<b>Se - 78</b>	17.8	64.0	60.5	5.5	115
<b>Ag - 107</b>	<0.003	34.1	34.6	1.6	85.2
<b>Cd - 111</b>	<0.01	40.7	40.5	0.5	102
<b>Tl - 205</b>	1.10	44.1	44.3	0.5	107
<b>Pb - 208</b>	<0.02	38.7	39.3	1.5	96.7

- ▶ FGD is one of the most difficult matrices for ICP-MS due to extremely high TDS, sulfates, alkali-alkaline earth metals
- ▶ Must overcome sample transport effects, signal drift, suppression/enhancement effects, and polyatomic interferences
- ▶ Discrete, fast sampling minimizes solids deposition on cones per sample while speeding up analysis time
- ▶ Collision/Reaction cell technology can effectively handle polyatomic interferences
- ▶ Combined instrument, accessories, and method conditions provide stable signal over many hours for accurate and precise analysis of FGD

