

# Trace Analysis of Compounds in Potable Water by On-Line SPE-LC/MS/MS

Joan Stevens, Ph.D.  
Sample Preparation Applications Scientist  
CSD  
NEMC 2012, Thursday August 9th

# Why the Interest in Emerging Contaminants (EC) and Pharmaceuticals and Personal Care Products (PPCPs) in Water?

- > 3,000 drugs and PPCPs
- Interest and fate of drugs, hormones
  - Hormones affecting fertility, changes in gender
- Limited knowledge of treatability for all contaminants
- High public visibility, concern



# Requirements for Emerging Contaminant Analysis Methods

- Selectivity/Specificity
- Sensitivity-MRL determination by LCMS/MS
- Target List
- QA/QC
- Ruggedness/matrix effects
- Cost

# Innovative LC-MS/MS Approach to EC Analysis

- Cost effective screen for as many analytes as possible
- Reporting levels- 5 to 10 ng/L (ppt)
- Known target analytes and identification from experts
- Accuracy and precision comparable to existing drinking water methods- 70-130% recovery
- Use multiple ion transitions (MRM)-misidentification

# Examples of EC Concerns

<b>Pharmaceuticals</b>	
Veterinary and human antibiotics	Trimethoprim, erythromycin, lincomycin, sulfamethaxazole
Analgesics, anti-inflammatory drugs	Codein, ibuprofene, acetaminophen, acetylsalicylic acid, diclofenac, fenoprofen
Psychiatric drugs	Diazepam
Lipid regulators	Bezafibrate, clofibrac acid, fenofibrac acid
β-blockers	Metoprolol, propanolol, timolol
X-ray contrasts	Iopromide, iopamidol, diatrizoate
<b>Steroids and hormones</b>	
Personal care products	Estradiol, estrone, estriol, diethylstilbestrol
Fragrances	Nitro, polycyclic and macrocyclic musks,
Sun-screen agents	Benzophenone, methylbenzylidene camphor
Insect repellents	N,N-diethyltoluamide
Antiseptics	Triclosan, Chlorophene
<b>Surfactants and surfactant metabolites</b>	
	Alkylphenol ethoxylates, 4-nonylphenol, 4-octylphenol, alkylphenol carboxylates
<b>Flame retardants</b>	
	Polybrominated diphenyl ethers (PBDEs),
	Tetrabromo bisphenol A, C <sub>10</sub> -C <sub>13</sub> chloroalkanes
	Tris (2-chloroethyl)phosphate
<b>Industrial additives and agents</b>	
	Chelating agents (EDTA), aromatic sulfonates,
<b>Gasoline additives</b>	
	Dialkyl ethers, Methyl- <i>t</i> -butyl ether (MTBE)

**Pesticides/Herbicides  
Sugar Substitutes**

# Sample Preparation Requirements: Basic Steps


- Collection device
- Collection of water sample 1000 mL
- Shipping/cost
- Analysis



# Sample Preparation

## Off-line SPE


- 1000 mL transported to the lab
- Trace enrich onto a cartridge or disk
- Elute with few mLs solvent, dry, reconstitute
- Amenable LC, % organic solvent
- 1 mL eluent loaded into autosampler vial
- 2-20  $\mu\text{L}$  injected, UHPLC to HPLC UV/MS detection
- End result 98-99% prepared sample is wasted



# Sample Preparation

## On-line SPE (trace enrichment-SPE)


- 100% of the prepared sample is loaded
- Volume can be <5 mL
- Combined with more sensitive detection (MS/MS)



- Logistics of getting sample to the lab is simplified



# Setup of the Agilent AQUA Online-SPE System



## SPE system:

- Quaternary pump
- G1329A autosampler with 900  $\mu$ l head
- 6 port 2 pos valve
- 12 port 6 pos valve
- 12 port stream selection valve (optional)

## RRLC system:

- Binary pump SL
- Well plate sampler SL (not used)
- Column department SL

# Method Implementation in Mass Hunter Acquisition

**Agilent MassHunter Workstation Data Acquisition**

Context: Acquisition | Layout: Default | Method: TG\_Online\_SPE\_EA\_w\_SPE\_ALS\_short\_inj3.m | Worklist: EA\_SPE\_26sept08-2.wkl

**Devices**

- ALS: 900  $\mu$ l, # 112
- BinPump-SL: 0.3 ml/min, A1 95%, B1 5%, 146.15 bar
- Quat Pump: 2 ml/min, A 100%, B 0%, C off, D off, 28.43 bar
- Column-SL: A 39.98 °C, B 40.01 °C
- Valve: G1158A, Current Pos: 2
- Valve: G1159A, Current Pos: 1
- MS QQQ: ESI(+), MS, Cal Off, CC Gas On

**Actuals**

Parameter	Value
Binary Pump: Ripple	3.15 %
QQQ: Rough Vac	1.97E+0 Torr
QQQ: High Vac	2.51E-5 Torr
QQQ: Turbo Speed	100.0 %
QQQ: MS1 Heater	100 °C
QQQ: MS2 Heater	100 °C
QQQ: Gas 1 temp	320 °C
QQQ: Vaporizer	29 °C
QQQ: Gas Flow	8.0 l/min
QQQ: Not Ready Text Short	
QQQ: Not Ready Text Long	
ALS: Temp	off

**Chromatogram (TIC)**: Binary Pump: Pressure (red), Quaternary Pump: Pressure (green). X-axis: min (-30 to 20). Y-axis: Abundance (0 to 1500).

**Mass Spectrum (QQQ Spectrum)**: MS [7: MRM (203.0), ESI (+), 17.36]. X-axis: m/z (amu) (110 to 170). Y-axis: Abundance (0 to 6). Peaks at m/z 104.0 and 175.0.

**Setup Panel**

- Flow: 2 ml/min
- Solvents: A: 100 % H<sub>2</sub>O, B: 0 % AcN, C: Off, D: Off
- Stop Time: No Limit, 15 min
- Post Time: Off, 0 min
- Pressure Limits: Min: 0 bar, Max: 200 bar

## Single 900 $\mu$ L injection:

DRAW def. amount from sample

VALVE mainpass

WAIT 0.7 min / 1.5 min\*\*

REMOTE startpulse

\*\* The delay volume of the online-SPE system is approximately 100  $\mu$ L. If the 900  $\mu$ L seat extension is used the wait time has to be extended to 1.5 min.



# On-line SPE method – Autosampler Program(s)

## Single 900 $\mu$ L injection:

DRAW def. amount from sample
------------------------------

VALVE mainpass
----------------

WAIT 0.7 min / 1.5 min**
--------------------------

REMOTE startpulse
-------------------

\*\* The delay volume of the system until the cartridge is approximately 100  $\mu$ L. If the 900  $\mu$ L seat extension is used the wait time has to be extended to 1.5 min (flush volume 3 mL).



## 2-fold 900 $\mu$ L injection using the multiple draw kit (900 $\mu$ L seat extension):

DRAW def. amount from sample
------------------------------

EJECT def. amount into seat
-----------------------------

DRAW def. amount from sample
------------------------------

VALVE mainpass
----------------

WAIT 1.5 min
--------------

REMOTE startpulse
-------------------



# Online-SPE Method – Loading Pump

Mobile phase:	A: 0.05% formic acid B: acetonitrile C: isopropanol			
SPE cartridge**:	PLRP-s available 2.1 x 12.5mm or 4.6 x 12.5mm			
Temperature:	ambient			
Flow (load):	2 mL/min			
Loading and equilibration:	0.00 min	0 % B	0 % C	2.0 ml/min
	0.10 min	0 % B	0 % C	0.1 ml/min
	12.50 min	0 % B	0 % C	0.1 ml/min
	12.60 min	0 % B	100 % C	2.0 ml/min
	12.70 min	0 % B	100 % C	2.0 ml/min
	14.00 min	0 % B	100 % C	2.0 ml/min
	14.10 min	100 % B	0 % C	2.0 ml/min
	15.50 min	100 % B	0 % C	2.0 ml/min
	15.60 min	0 % B	0 % C	2.0 ml/min



\*\* SPE cartridges fit in the Agilent Guard Column Hardware Kit (P.N. 820888-901)  
All connections between the analytical system and the SPE hardware have been made by using  
0.18 mm ID PEEK tubing (P.N. 5042-6462) and finger-tight PEEK fittings (P.N. 5063-6591).

# Online-SPE – Timetable for Loading and Analysis



Equilibration of SPE cartridges, SPE and separation of pesticides is overlapped and results in a total runtime for the method of 20 minutes.

# Example 1: High Sensitivity Analysis of Pesticides

□ Dqdd }hg#khuelflgh#frp srxqgv



Atrazine

Carbendazim

Carbetamid

Chloridazon

Chlorotoluron

Chloroxuron

Desethylatrazine

Desisopropylatrazine

Desmetryn

Diflubenzuron

Diuron

Fenuron

Irgarol 1051

Isoproturon

Linuron

Metamitron

Methabenzthiazuron

Metoxuron

Metsulfuron methyl

Monolinuron

Monuron

Neburon

Prometryn

Propazin

Simazine

Terbutryn

Terbutylazine

Trietazin

# Example 1: High Sensitivity Analysis of Pesticides

□ V | vhp #Erqilj xudwlrq



AQUA Online-SPE system consisting of

- Quarternary pump
- G1329A autosampler with 900 µl head
- 6 port 2 pos valve
- 12 port 6 pos valve



UHPLC system consisting of



- 1290 Infinity binary pump
- 1290 Infinity autosampler (not used)
- 1290 Infinity thermostatted column compartment



G6460AA Triple Quadrupole Mass Spectrometer

# Example 1: High Sensitivity Analysis of Pesticides

□ Inversurkung





- Quantifier: 207.1  $\rightarrow$  72.0
- Qualifier: 207.1  $\rightarrow$  46.1
- Peak width: 16.1 s
- Data points: 18
- Linear range: < 1 ng/L to 1  $\mu$ g/L






# Example 1: High Sensitivity Analysis of Pesticides

☐ Ivrsurwurg +V2Q #fdo #shdn Cw Cshdn,



# Example 1: High Sensitivity Analysis of Pesticides



☐ F | d q d } l q h



- Quantifier: 241.1 → 214.0
- Qualifier: 241.1 → 68.0
- Peak width: 16.0 s
- Data points: 18
- Linear range: 1 ng/L to 5 µg/L




# Example 1: High Sensitivity Analysis of Pesticides

□ F | d q d } l q h + V 2 Q # F d o F # s h d n C w C s h d n ,



# Example 1: High Sensitivity Analysis of Pesticides



Ehqw d } r q h



- Quantifier: 239.1  $\rightarrow$  132.0
- Qualifier: 239.1  $\rightarrow$  197.0
- Peak width: 19.8 s
- Data points: 23
- Linear range:  $\ll$  1 ng/L to 0.1  $\mu$ g/L


# Example 1: High Sensitivity Analysis of Pesticides

□ E h q w d } r q h + V 2 Q # f d d # s h d n ( w r C s h d n ,



# Example 1: High sensitivity analysis of pesticides

53Cirog#qnfwrq#r#Fdb7#83#bj20#vdp h#Edwulgjh,



Compound	Average (µg/L)	SD	RSD
2OH-Atrazine	0.059	0.002	3.50
Atrazine-desisopropyl	0.051	0.001	2.60
2,6-Dichlorobenzamide	0.047	0.003	6.80
Metamitron	0.053	0.001	2.64
Atrazine-desethyl	0.051	0.002	3.43
Dimethoate	0.055	0.002	2.94
Hexazinone	0.060	0.002	2.56
Simazine	0.052	0.001	2.84
Cyanazine	0.043	0.001	2.79
Bentazone	0.057	0.002	3.32
Atrazine	0.053	0.001	2.34

Compound	Average (µg/L)	SD	RSD
Isoproturon	0.053	0.001	2.56
2,4-D	0.055	0.002	3.71
MCPA	0.046	0.003	5.97
DNOC	0.044	0.003	6.31
Mechlorprop	0.050	0.002	3.98
Dichlorprop	0.047	0.003	5.62
Terbutylazin	0.053	0.001	2.45
Dinoseb	0.039	0.002	4.55
Pendimethalin	0.040	0.002	4.54

## Example 2: Impact of Solvent and Sample Composition


□ Uhtx.luhp hqww

- In addition to method for the neutral herbicides an on-line SPE method for the following 5 compounds (LOQ < 10 ng/L):
  - Flufenacet
  - Iodsulfuron-methyl
  - Isoxaflutole
  - Mesosulfuron-methyl
  - Mesotrione
- Sample solution contains 80 mg/L Na<sub>2</sub>SO<sub>4</sub> for sample conservation
- Compounds belong to group of acidic herbicides, Mesotrione shows broad peak in chromatography

# Example 2: Impact of Solvent and Sample Composition

☐ K.SOF#Ekrp dwrj udp


- Calibration standard 5 (110 to 197 ng/L); injection volume 900  $\mu$ L





## Example 2: Impact of Solvent and Sample Composition


- Flow rate of the loading pump had no influence on loading efficiency and therefore a flow of 2 mL/min was chosen to reduce loading time.
- Compared to the use of demineralized water as the loading solvent loading with 0.1% formic acid increased the response significantly.



- For injection volumes > 900  $\mu$ L the addition of 0.1% formic acid to the sample solution is strongly recommended.

## Example 2: Impact of Solvent and Sample Composition

- The effect of the addition of acetonitrile and  $\text{Na}_2\text{SO}_4$  (80 mg/L) to the sample solution on the on-line SPE recovery was tested




- Acetonitrile addition was detrimental, the addition of  $\text{Na}_2\text{SO}_4$  increased the mesosulfurone-methyl response. It is recommended to prepare samples and calibration standards identical.

# Example 3: Robust Method for Glyphosate and AMPA

## □ Edfnjurxqg

- Glyphosate is a global herbicide which is widely used in agriculture and urban landscape management
- In the environment glyphosate is metabolized to its metabolite aminomethyl phosphonic acid (AMPA)
- Both compounds are extremely polar due to their bipolar structure
- For the analysis of both compounds a derivatization is widely accepted



## Example 3: Robust Method for Glyphosate and AMPA


### □ Edfnjurxqg

- Derivatization works at basic pH-values with high excess of FMOC and both reduces the lifetime of the chromatographic column and the robustness of the method
- Sensitivity of direct injection is just enough to detect 100 ng/L
- Online-SPE allows for clean-up and enrichment to comply with screening
- Experimental:
  - 1 mL of sample is fortified with 100 ng ISTD (Glyphosate 1,2 <sup>13</sup>C-<sup>15</sup>N, AMPA <sup>13</sup>C)
  - Borate buffer (pH 10) and FMOC (5 mg/mL) is added and reaction takes place over night (> 4 h)
  - Acetic acid (58%) is added to neutralize sample
  - FMOC derivates are stable for > 48 hours


# Example 3: Robust Method for Glyphosate and AMPA

DP SD


AMPA - 7 Levels, 7 Levels Used, 20 Points, 20 Points Used, 0 QCs



+ MRM (334.0 -> 179.0) Std\_10-r001.d



334.0 -> 179.0 , 334.0 -> 112.0



- Quantifier: 334.0  $\rightarrow$  179.0
- Qualifier: 334.0  $\rightarrow$  112.0
- Peak width: 15.4 s
- Data points: 15
- Linear range: 10 ng/L to 1.0  $\mu$ g/L

## Example 3: Robust Method for Glyphosate

### □ UHVXOW

- Use of internal standards needed to correct for derivatization yields in different water samples
- Online-SPE increases the sensitivity to the relevant concentration range (10 to 25 ng/L) for environmental samples
- Robustness of the method is substantially increased – run already > 1000 samples on the same chromatographic column with just 2 cartridges
- AMPA has been found in concentrations up to 2000 ng/L, Glyphosate just positive if AMPA concentrations are extremely high

## Review: On-line versus Off-Line Trace Enrichment (SPE)

On-Line	Off-Line
Increased sensitivity-inject entire sample-analyzed	Only a fraction of concentrated sample injected 2-100 uL from 1 mL extract
Small Sample Volumes (mL)	Large sample volume 100mL-L
Increased sample throughput	Extraction and analysis are separate
Reduce error-sample handling	Increase chance of error-manual
No evaporation step	Concentration step needed
No reconstitution	Reconstitution
Integrated System	Standalone Instruments

# Summary and Conclusions

- The on-line SPE combined with a LC-QQQ system allows a relatively simple, fast and reliable determination of herbicides in the low ng/L range in filtered water samples.
- The whole system is fully controlled with the MassHunter acquisition software.
- Adding a more sensitive QQQ system allows for even better sensitivity or compounds which are less weakly ionizable.
- Good recovery values and reproducibilities can be achieved even in complex samples and for very polar compounds.
- Online SPE not only increases the sensitivity but adds robustness to the method (e.g. for glyphosate)
- The chemistry remains the same – adding acid to the loading solvent increases the recovery for acidic compounds, adding salt to the sample solution increases the adsorption coefficient for some compounds.
- The use of online SPE with scanning instruments allows for unknown screening in complex environmental samples.





**I'M GOING IN  
FOR A SAMPLE**

***Thank You!  
Any Questions?***