# Concerns Regarding 24-hr Sampling for Formaldehyde, Acetaldehyde, and Acrolein with 2,4-dinitrophenylhydrzine (DNPH)-coated Solid Sorbents

#### Jason Herrington, Ph.D. Air Innovations Chemist





Innovative Chromatography Products www

#### The Five Ws

- Who
- What
- When
- Where
- Why







Innovative Chromatography Products www

# WHAT... are DNPH-coated solid sorbents?

- Silica-Gel packed cartridge coated with acidified 2,4-dinitrophenyldhydrzine (DNPH)
- Air is drawn over the cartridge actively or passively
  - Carbonyls (aldehydes and ketones) are collected by derivatizing with DNPH on cartridge
- Sample is extracted with ACN and analyzed with HPLC-UV technique









Innovative Chromatography Products

### WHERE... are the used?

#### • The list goes on ad infinitum...











Innovative Chromatography Products

## WHO... wrote this method?



- Several national and international standardized methods
  - U.S. EPA Method TO-11A
  - CARB 1004
  - ASTM 5197
  - Deutsche Forschungsgemeinschaft Methoden Nr.2

#### "Gold Standard"





### WHAT... concerns?

Method has been around for 30 years

Extensively developed and evaluated

- What concerns... what could have been missed?
  - Wrong question!



"YOU REALLY SHOULD GIVE ME A LITTLE CREDIT, NOT MANY MEN COULD COPE WITH FAILURE THE WAY I HAVE."

What has not been accurately reflected in the standardized methods?





## WHAT... compounds?

- DNPA interference with Formaldehyde
  - Known human carcinogen



ED PERSONNEL ONLY

- Low collection efficiencies for 24-hour sampling Acetaldehyde
  - Probable human carcinogen
- Plagued by several issues for Acrolein
  - Highly actuely toxic







Innovative Chromatography Products

#### WHAT... is DNPA?

- DNPH reacts with NO<sub>2</sub> to form 2,4-Dinitrophenylazide (DNPA)
  - DNPA coelutes with the formaldehyde-DNPH derivative





Innovative Chromatography Products www.



Adopted from Potter and Karst, 1996





Innovative Chromatography Products

## WHAT... is the solution for DNPA?

- Altering/optimizing the HPLC gradient can overcome the DNPA interference.
  - Several studies have capitalized on the formation of DNPA for the sampling of NO<sub>2</sub>, by altering/optimizing their HPLC gradient

Diode array detector set at 300 nm for DNPA can be utilized





Innovative Chromatography Products www.restek.com

# WHY... am I talking about DNPA then?

Remember OZONE???

The use of KI scrubbers (as recommended in the methods) to limit  $O_3$  interferences promotes  $NO_2$  by the oxidation of NO.

The standardized methods do not reflect any of the aforementioned information

– TO-11A 1999 revision overlooked this issue





Innovative Chromatography Products www

### WHAT... is up with acetaldehyde

- Through an extensive literature search, we were only able to find three studies during which carbonyls other than formaldehyde were evaluated on DNPH-coated solid sorbents for long-term sampling (i.e., 24 hour or greater).
- Lazarus (1999) reported low acetaldehyde collection efficiencies (CE); and Grosjean (1991), and Grosjean and Grosjean (1995) evaluated breakthrough of the collection media, which does not necessarily reflect CE.





Innovative Chromatography Products www

| Experimental condition | Carbonyl     | SUPELCO                 | WATERS              | XPOSURE             | HOUSE               |
|------------------------|--------------|-------------------------|---------------------|---------------------|---------------------|
| 3 hours at 30% RH      | Formaldehyde | $89 \pm 10^{\circ}$ (3) |                     |                     |                     |
|                        | Acetaldehyde | 93 ± 8° (3)             |                     |                     |                     |
| 24 hours at 30% RH     | Formaldehyde | 83 ± 4<br>(3)           | 87 ± 11<br>(3)      | 111 ± 4<br>(3)      | $104 \pm 25$ (3)    |
|                        | Acetaldehyde | 39 ± 7<br>(3)           | 43 ± 3<br>(3)       | 62 ± 7<br>(3)       | 1 ± 2<br>(3)        |
| 48 hours at 30% RH     | Formaldehyde | 89 ± 8<br>(3)           | 93 ± 4<br>(3)       | 105 ± 19<br>(3)     | 14 ± 8<br>(3)       |
|                        | Acetaldehyde | $51 \pm 22$ (3)         | 43 ± 2<br>(3)       | 40 ± 11<br>(3)      | 0<br>(3)            |
| 24 hours at 60% RH     | Formaldehyde | 101 ± 8<br>(3)          | $101 \pm 13$<br>(3) | $121 \pm 32$<br>(3) | $133 \pm 27$<br>(3) |
|                        | Acetaldehyde | 27 ± 4<br>(3)           | 29 ± 2<br>(3)       | 30 ± 2<br>(3)       | 9 ± 2<br>(3)        |

Herrington, J.; Fan, Z.; Lioy, P. J.; Zhang, J. Low acetaldehyde collection efficiencies for 24-hour sampling with 2,4dinitrophenylhydrazine (DNPH)-coated solid sorbents. Environ. Sci. Technol. 2007, 41, 580-585.







Herrington, J.; Fan, Z.; Lioy, P. J.; Zhang, J. Low acetaldehyde collection efficiencies for 24-hour sampling with 2,4dinitrophenylhydrazine (DNPH)-coated solid sorbents. Environ. Sci. Technol. 2007, 41, 580-585.





# WHAT... is the solution for acetaldehyde?

To date... there is none for DNPH-based methods

#### Again... the standardized methods do not reflect any of the aforementioned information





Innovative Chromatography Products www.re



- •L. Benning and A. Wahner, J. Atmos. Chem., 1998, 31, 105–117.
- •E. Goelen, M. Lamrechts and F. Geyskens, Analyst, 1997, 122, 411–419.
- •S. B. Tejada, J. Environ. Anal. Chem., 1986, 26, 167–185.
- •K. Olson and S. J. Swarin, J. Chromatogr., 1985, 333, 337–347.
- •M. Possanzini and V. DiPalo, Chromatographia., 1995, 40,134–138.
- •C. H. Risner and P. Martin, J. Chromatogr. Sci., 1994, 32, 76–82.
- •C. H. Risner and P. Martin, J. Chromatogr. Sci., 1994, 32, 76–82.
- •A. Sakuragawa, T. Yoneno, K. Inoue and T. Okutani, J. Chromatogr., 1999, 844A, 403–408.
- •R. Schulte-Ladbeck, R. Lindahl, J. O. Levin and U. Karst, J. Environ. Monit., 2001, 3, 306–310.





#### WHAT... is the solution

- In 2010 acrolein was finally removed from the list of TO-11A target analytes
- I hear rumors of DNPH methods, but yet to see any methods published in the peer reviewed literature.
- Other derivatizing agents:
  - DNSH
  - PFBHA





Innovative Chromatography Products w

# WHY... am I talking about old news?

 Current methods do not reflect the aforementioned information

Therefore, end users are completely un-aware of the aforementioned

To date... no single DNPH-based method has addressed all of the issues outlined here





Innovative Chromatography Products www.re

#### WHEN... did all this take place?

- First discovered:

   Acrolein 1985
   2010 (25 years later)

   Formaldehyde 1996
   ???? (16 years so far)
  - Acetaldehyde 2007 ???? (5 years so far)





Innovative Chromatography Products w

# WHERE... do we go from here?

- U.S. EPA (and other standardizing entities) need to publish appropriate amendments to the methods NOW
  - End users need to be made aware
  - Laboratories need to develop methods for formaldehyde, acetaldehyde, and acrolein... at a minimum
    - No more narrow-sightedness on just acrolein
- Need to rigorously test methods
  - Use a dynamic gas generation system
  - Evaluate temperature, relative humidity, time, etc... "Real World" variables
  - No more assumptions!





Innovative Chromatography Products www.restek.com

#### WHAT... more questions?

- ChromaBLOGraphy
  - http://blog.restek.com/
- Herrington and Hays, 2012. Concerns regarding 24-h sampling for formaldehyde, acetaldehyde, and acrolein using 2,4-dinitrophenylhydrazine (DNPH)-coated solid sorbents. Atmospheric Environment, 55, 179-18.





Innovative Chromatography Products www.restek.com