

Presented by: Stephen T. Zeiner, CEAC Co-author: Kristina Harsh August 2012 National Environmental Monitoring Conference



Summary

- Why evaluate contamination?
- Sources of Contamination
- Validation Guidelines
- Examples of Contamination Assessments
- Comparability Issues
- Path Forward





- Method sensitivity/selectivity
 - False positives
 - Interferences
- Contamination/qualification impacts:
 - Notice of Violation
 - Additional discharge sampling
 - Additional remediation
 - Not enough remediation
 - Improper risk assessment





- Collection equipment
- Bottleware
- Preservatives
- Preparation equipment
- Storage
- Instrument carryover
- Reagents (solvents/spikes)



October 1999 National Organic Validation Guidelines

Blank Result	Sample Result	Action for Samples
Detect	>CRQL but <5×blank result*	Qualify results as "U"
Detect	<crql and<br=""><5×blank result*</crql>	Report CRQL value with a "U"
Detect	>5×blank result*	No qualification
Gross contamination	Detects	Qualify results as unusable "R"

*10×the CRQL for common laboratory contaminants



June 2008 National Organic Validation Guidelines

Blank Result	Sample Result	Action for Samples
<crql< td=""><td><crql*< td=""><td>Report CRQL value with a "U"</td></crql*<></td></crql<>	<crql*< td=""><td>Report CRQL value with a "U"</td></crql*<>	Report CRQL value with a "U"
<crql< td=""><td>≥CRQL*</td><td>Use professional judgment</td></crql<>	≥CRQL*	Use professional judgment
>CRQL	<crql*< td=""><td>Report CRQL value with a "U"</td></crql*<>	Report CRQL value with a "U"
>CRQL	≥CRQL* and <blank concentration<="" td=""><td>Report the blank concentration for the sample with a "U" or qualify data as unusable "R"</td></blank>	Report the blank concentration for the sample with a "U" or qualify data as unusable "R"

*2×CRQL for methylene chloride, 2-butanone, and acetone or 5×CRQL for bis(2-ethylhexyl)phthalate



June 2008 National Organic Validation Guidelines

Blank Result	Sample Result	Action for Samples
>CRQL	≥CRQL* and ≥blank result	Use professional judgment
=CRQL	<crql*< td=""><td>Report CRQL with a "U"</td></crql*<>	Report CRQL with a "U"
=CRQL	≥CRQL*	Use professional judgment
Gross contamination	Detect	Qualify results as unusable "R"

*2×CRQL for methylene chloride, 2-butanone, and acetone or 5×CRQL for bis(2-ethylhexyl)phthalate



February 1994 National Inorganic Validation Guidelines

Blank Result	Sample Result	Action for Samples
Detect	<5×blank result	Qualify results as "U"
Detect	>5×blank result	No qualification
Gross contamination	Detects	Qualify results as unusable "R"



October 2004 National Inorganic Validation Guidelines

Calibration Blanks

Blank Result	Sample Result	Action for Samples
≤CRQL	≤CRQL	Report CRQL with a "U"
≤CRQL	>CRQL	Use professional judgment
>CRQL	≤CRQL	Report CRQL with a "U"
>CRQL	>CRQL but <blank result</blank 	Report at level of blank result with a "U" or qualify data as unusable "R"
>CRQL	>blank result	Use professional judgment



October 2004 National Inorganic Validation Guidelines

Preparation Blanks

Blank Result	Sample Result	Action for Samples
≤CRQL	≤CRQL	Report CRQL with a "U"
≤CRQL	>CRQL	Use professional judgment
>CRQL	≤CRQL	Report CRQL with a "U"
>CRQL	>CRQL but <10×blank result	Qualify results as unusable "R" or estimated "J"
>CRQL	>10×blank result	No qualification





Blank Result	Sample Result	Action for Samples
Detect	Not detected	No qualification
<crql< td=""><td><ql <5×blank<br="" and="">result*</ql></td><td>Report CRQL value with a "U"</td></crql<>	<ql <5×blank<br="" and="">result*</ql>	Report CRQL value with a "U"
<crql< td=""><td>>QL but <5×blank result*</td><td>Qualify results as "U"</td></crql<>	>QL but <5×blank result*	Qualify results as "U"
<crql< td=""><td><ql and="">5×blank result*</ql></td><td>No further action</td></crql<>	<ql and="">5×blank result*</ql>	No further action
<crql< td=""><td>>QL and >5×blank result*</td><td>No qualification</td></crql<>	>QL and >5×blank result*	No qualification

* Common contaminants would be 10×blank result.



Region 1 Organic Validation Guidelines

Blank Result	Sample Result	Action for Samples
≥CRQL	>QL but ≤5×blank result*	Report value with a "U"
≥CRQL	<ql and="" ≤5×blank<br="">result*</ql>	Qualify results as "U"
≥CRQL	≥ QL and >5×blank result*	No further action
≥CRQL	<ql but="">5×blank result*</ql>	No qualification
Gross contamination (> 10 × CRQL)	Detects	Qualify results as unusable "R"





Blank Result	Sample Result	Action for Samples
≥MDL	≥QL and >5×blank result	No qualification
≥MDL	<ql and="">5×blank result</ql>	No further action (report estimated sample result)
≥MDL	≥QL but ≤5×blank result	Qualify results as "U"
≥MDL	< QL and \leq 5 × blank contaminate	Report CRQL value with a "U"
≥MDL	Non-detect (U)	No qualification



Region 2 Organic Validation Guidelines

Blank Result	Sample Result	Action for Samples
Detects	Not detected	No qualification
<crql< td=""><td><crql< td=""><td>Report CRQL value with a "U"</td></crql<></td></crql<>	<crql< td=""><td>Report CRQL value with a "U"</td></crql<>	Report CRQL value with a "U"
<crql< td=""><td>≥CRQL</td><td>Use professional judgment</td></crql<>	≥CRQL	Use professional judgment
>CRQL	≥CRQL and <blank result<="" td=""><td>Report CRQL value with a "U"</td></blank>	Report CRQL value with a "U"
>CRQL	≥CRQL and ≥blank result	Use professional judgment
=CRQL	<crql< td=""><td>Report CRQL value with a "U"</td></crql<>	Report CRQL value with a "U"
=CRQL	≥CRQL	Use professional judgment
Gross contamination	Detects	Qualify results as unusable "R"





Instrument Blanks

Blank Result	Sample Result	Action for Samples
≤CRQL	≤CRQL	Report CRQL value with a "U"
>CRQL	>CRQL but <blank result</blank 	Qualify results as unusable "R"
>CRQL	>CRQL and ≤10×blank result	Qualify results as estimated "J"
>CRQL	≤CRQL	Report CRQL value with a "U"
<(-CRQL)	≥CRQL but <10×CRQL	Qualify results as estimated "J"





Preparation Blanks

Blank Result	Sample Result	Action for Samples
≤CRQL	≤CRQL	Report CRQL value with a "U"
>CRQL	>Blank result but <10×blank result	Flag results as estimated "J"
>CRQL	>CRQL but <blank p="" result<=""></blank>	Reject "R" all sample results
>CRQL	Prep blank < field blank	No qualification
>CRQL	≤CRQL	Report CRQL value with a "U"



Region 3 Organic Validation Guidelines

Blank Result	Sample Result	Action for Samples
Detected	<5×blank result*	Qualify results as "B"
Detected	Not-detected or >5×blank result*	No qualification
Gross contamination	Detects	Qualify results as unusable "R"

* Common contaminants would be 10×blank result.





Blank Result	Sample Result	Action for Samples
Detected	<5×blank result	Qualify results as "B"
Detected	Not-detected or >5×blank result	No qualification

The blank result is the highest concentration of all associated blanks.





- Method blank has 8.2 µg/L reported
- Quantitation limit is 200 µg/L
- Sample has 85 µg/L reported
- NFG 1999 qualification: None
- NFG 2008 qualification: 200 U
- Region 1 qualification: None
- Region 2 qualification: 200 U
- Region 3 qualification: None





- Method blank has 27.0 µg/L reported
- Quantitation limit is 25.0 µg/L
- Sample has 190 µg/L reported
- NFG 1999 qualification: 190 U
- NFG 2008 qualification: Use prof. judgment
- Region 1 qualification: 190 U
- Region 2 qualification: Use prof. judgment
- Region 3 qualification: 190 B





- Method blank has 0.013 µg/L reported
- Quantitation limit is 0.020 µg/L
- Sample has 0.023 µg/L reported
- NFG 1999 qualification: 0.023 U
- NFG 2008 qualification: None
- Region 1 qualification: 0.023 U
- Region 2 qualification: Use prof. judgment
- Region 3 qualification: 0.023 B





- Method blank has 1.5 µg/L reported
- Quantitation limit is 0.30 µg/L
- Sample has 6.0 µg/L reported
- NFG 1994 qualification: 6.0 U
- NFG 2004 qualification: 6.0 J or 6.0 R
- Region 1 qualification: 6.0 U
- Region 2 qualification: 6.0 J
- Region 3 qualification: 6.0 B





- Method blank has 8.5 µg/L reported
- Quantitation limit is 5.0 µg/L
- Sample has 7.0 µg/L reported
- NFG 1994 qualification: 7.0 U
- NFG 2004 qualification: 7.0 U or 7.0 R
- Region 1 qualification: 7.0 U
- Region 2 qualification: 7.0 R
- Region 3 qualification: 7.0 B



Comparability Issues

- Use of different criteria for contamination evaluation is a kin to changing data quality objectives in the middle of a project.
- Differences in approaches to contamination evaluation can result in problems with:
 - Continuity of data over time
 - Comparison of data across regions/states
 - Comparison of validation between firms
 - Data use for risk assessment



Path Forward

- Need to address:
 - Difference in Guidelines
 - Professional Judgment
- Ways to address:
 - Quality Assurance Project Plan
 - Standard Operating Procedures
 - Reports



References

- Region I, EPA-New England Data Validation Functional Guidelines for Evaluating Environmental Analyses, July 1996, Revision December 1996
- USEPA Hazardous Waste Support Branch Validating Volatile Organic Compounds by Gas Chromatography/Mass Spectrometry SW-846 method 8260B, SOP #HW-24, Revision #2, August 2008
- USEPA Hazardous Waste Support Branch Validating Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry SW-846 Method 8270D, SOP #HW-22, Revision #4, August 2008



References

- USEPA Hazardous Waste Support Branch Validation of Metals for the Contract Laboratory Program (CLP) based on SOW ILM05.3, SOP #HW-2, Revision #13, September 2006
- Region III Modifications to National Functional Guidelines for Organic Data Review Multi-Media, Multi-Concentration (OLMO1.0-OLMO1.9), September 1994
- Region III Modifications to the National Functional Guidelines for Evaluating Inorganics Analyses, April 1993
- USEPA Contract Laboratory Program National Functional Guidelines for Superfund Organic Methods Data Review, October 1999 and June 2008
- USEPA Contract Laboratory Program National Functional Guidelines for Inorganic Data Review, February 1994 and October 2004







Environmental Standards, Inc. "Setting the Standards for Innovative Environmental Solutions"

Headquarters 1140 Valley Forge Road | PO Box 810 | Valley Forge, PA 19482 | 610.935.5577
Virginia 1208 East Market Street | Charlottesville, VA 22902 | 434.293.4039
Tennessee 1013 Brentwood Way | Kingston, TN 37763 | 865.376.7590
Texas 11200 Richmond Avenue, Suite 350 | Houston, TX 77082 | 281.752.9782
Web www.envstd.com | E-mail solutions@envstd.com

